


**Revision Plate 2.7mm-3.5mm
Surgical Technique**

Astrolabe recognizes that proper surgical procedures and techniques are responsibilities of medical professionals.

The following guidelines are provided for information purposes only. Each surgeon must evaluate the appropriateness of the procedures based on their medical training, experience and condition of the patient. Before using the system, the surgeon must consult the operating instructions for additional warnings, precautions, indications, contraindications and adverse effects.

Revision Plate 2.7mm / 3.5mm



Revision Plate 2.7mm / 3.5mm



2.7 mm
Cortical/ Locking

3.5mm
Cortical/ Locking



Locking Plate, Revision, Left



Locking Plate, Revision, Right

General Features

- The anatomical shape of the revision plate allows its placement without excessive flexion or even flexion due to its correspondence with the first ray with the distal phalange of the big toe and the development of strength recovered from the FH tendon on the IP-joint.
- Distal fixation with 2.7 mm Locking/ Cortical Screws and proximal fixation with 3.5 mm Locking/ Cortical Screws
- 2.0mm thickness

Indications:

- Keller-Brandes Revisions;
- MTP 1 lengthening arthrodesis;
- Revisions after MTP 1 prosthesis;
- MT 1 lengthening osteotomy.

- After choosing the appropriate plate, provisionally fix it in place with a Kirschner Wire.
- If it is necessary, the plate can be moulded (item 12) to better fit patient's anatomy.



- **Note:** Avoid excessive modeling of the plate as this can compromise its locking mechanisms. When using plate bender (item 12), holes adjacent to the bender can lose the ability to lock. If this occurs, a Cortical Screw must be used.

● 3.5 mm Locking Screws Positioning

Properly place the Threaded Graduated Drill Guide (item 07) by turning it into the threaded hole, then proceed with drilling using the Drill Bit (item 02/04).



- The reading of the screw measurement can be made directly on the Threaded Graduated Drill Guide (item 07) or using the Depth Gauge (item 10).

- Use the Handle (item 05) and Shaft Screwdriver (item 09) to position the screws.



Attention:

To ensure the perfect fit and low profile of the Locking Screws on the plate, the Threaded Graduated Drill Guide (item 07) must be used.

If there is any damage to the threaded hole of the plate, the Locking Screws must be replaced by Cortical Screws.

● 2.7 mm Locking Screws Positioning

Properly place the Threaded Graduated Drill Guide (item 06) by turning it into the threaded hole, then proceed with drilling using the Drill Bit (item 01/03).



- The reading of the screw measurement can be made directly on the Threaded Graduated Drill Guide (item 06) or using the Depth Gauge (item 10).

- Use the Handle (item 05) and Shaft Screwdriver (item 08) to position the screws.




Attention:

To ensure the perfect fit and low profile of the Locking Screws on the plate, the Threaded Graduated Drill Guide (item 06) must be used.

If there is any damage to the threaded hole of the plate, the Locking Screws must be replaced by Cortical Screws.

- The placement of the screws procedure is repeated as many times as necessary, for optimal fixation of the plate.
- Check the final position of the screws through the image intensifier.
- Proceed with x-ray to check if final position is according to initial intention.




01

Drill Bit, Ø2.0 x 120 mm, Stop 50 mm,
AO Coupling, Blue Code
Cod.:09.01.03.20020



02

Drill Bit, Ø2.5 x 120 mm, Stop 50 mm,
AO Coupling, Yellow Code
Cod.:09.01.03.25020




03

Drill Bit, Ø2.0 x 125 mm, Stop 50 mm,
Stryker Coupling, Barrel Ø4.5 mm,
Blue Code
Cod.:09.01.07.20021



04

Drill Bit, Ø2.5 x 125 mm, Stop 50 mm,
Stryker Coupling, Barrel Ø4.5 mm,
Yellow Code
Cod.: 09.01.07.25021



05

Handle, Cannulated,
AO Coupling, 120 mm – Black
Cod.: 09.04.04.12050



06

Graduated Drill Guide
Ø2.0 x 40 mm, Threaded, Blue Code
Cod.:09.05.14.04020



07

Graduated Drill Guide
Ø2.5 x 40 mm, Threaded, Yellow Code
Cod.: 09.05.14.04025



08

Shaft Screwdriver, Torx-8
90 mm, AO Coupling, Blue Code
Cod.: 09.07.04.08091



09

Shaft Screwdriver, Torx-10
90 mm, AO Coupling
Yellow Code
Cod.:09.07.04.10092



10

Depth Gauge, 60 mm
Cod: 09.08.01.00060



11

Plate and Screw Holding
Forceps, Angled, 150 mm
Cod.09.10.06.00150



12

Bender
for 2.7/3.5 mm System Plates
Cod.09.13.00.02735





ASTROLABE
life and mobility

2.7mm- 3.5mm Revision Plate Surgical Technique

Rua José Gomes Ferreira nº 2 - Armazém 1
2660-517 São Julião do Tojal, Loures, Portugal
Tlf.: (+351) 219 672 298 | info@astrolabe-medical.com
www.astrolabe-medical.com