



Astrolabe

LIFE AND MOBILITY

Box Components

Hallux Valgus 2,7 & 3,5


Box Components Kit

09.50.03.00066


Component Cases


09.50.04.00066

09.50.04.00067



Hallux Valgus
2,7 & 3,5





Implant Case Set 1


Implant Case Set 2

N#	REF	DESC	QTY
1	15.04.03.27012	Cortical Screw, Ti, Ø2.7 x 12mm TX, S/I Extrem	5
	15.04.03.27014	Cortical Screw, Ti, Ø2.7 x 14mm TX, S/I Extrem	5
	15.04.03.27016	Cortical Screw, Ti, Ø2.7 x 16mm TX, S/I Extrem	5
	15.04.03.27018	Cortical Screw, Ti, Ø2.7 x 18mm TX, S/I Extrem	5
	15.04.03.27020	Cortical Screw, Ti, Ø2.7 x 20mm TX, S/I Extrem	5
	15.04.03.27022	Cortical Screw, Ti, Ø2.7 x 22mm TX, S/I Extrem	5
	15.04.03.27024	Cortical Screw, Ti, Ø2.7 x 24mm TX, S/I Extrem	5
2	15.04.03.27026	Cortical Screw, Ti, Ø2.7 x 26mm TX, S/I Extrem	5
	15.04.03.35012	Cortical Screw, Ti, Ø3.5 x 12mm TX, S/I Extrem	5
	15.04.03.35014	Cortical Screw, Ti, Ø3.5 x 14mm TX, S/I Extrem	5
	15.04.03.35016	Cortical Screw, Ti, Ø3.5 x 16mm TX, S/I Extrem	5
	15.04.03.35018	Cortical Screw, Ti, Ø3.5 x 18mm TX, S/I Extrem	5
	15.04.03.35020	Cortical Screw, Ti, Ø3.5 x 20mm TX, S/I Extrem	5
	15.04.03.35022	Cortical Screw, Ti, Ø3.5 x 22mm TX, S/I Extrem	5
	15.04.03.35024	Cortical Screw, Ti, Ø3.5 x 24mm TX, S/I Extrem	5
	15.04.03.35026	Cortical Screw, Ti, Ø3.5 x 26mm TX, S/I Extrem	5



Implant Cases Screw Sets

Implant Case Set 1



1


3

5

...

4

Implant Case Set 2





2

3


6











...


4

N#	REF	DESC	QTY
	15.14.03.27012	Locking Screw, Ti, Ø2.7 x 12mm TX, S/I Extrem	10
	15.14.03.27014	Locking Screw, Ti, Ø2.7 x 14mm TX, S/I Extrem	10
	15.14.03.27016	Locking Screw, Ti, Ø2.7 x 16mm TX, S/I Extrem	10
	15.14.03.27018	Locking Screw, Ti, Ø2.7 x 18mm TX, S/I Extrem	10
	15.14.03.27020	Locking Screw, Ti, Ø2.7 x 20mm TX, S/I Extrem	10
	15.14.03.27022	Locking Screw, Ti, Ø2.7 x 22mm TX, S/I Extrem	10
	15.14.03.27024	Locking Screw, Ti, Ø2.7 x 24mm TX, S/I Extrem	10
	15.14.03.27026	Locking Screw, Ti, Ø2.7 x 26mm TX, S/I Extrem	10
	15.14.03.35012	Locking Screw, Ti, Ø3.5 x 12mm TX, S/I Extrem	10
	15.14.03.35014	Locking Screw, Ti, Ø3.5 x 14mm TX, S/I Extrem	10
	15.14.03.35016	Locking Screw, Ti, Ø3.5 x 16mm TX, S/I Extrem	10
	15.14.03.35018	Locking Screw, Ti, Ø3.5 x 18mm TX, S/I Extrem	10
	15.14.03.35020	Locking Screw, Ti, Ø3.5 x 20mm TX, S/I Extrem	10
	15.14.03.35022	Locking Screw, Ti, Ø3.5 x 22mm TX, S/I Extrem	10
	15.14.03.35024	Locking Screw, Ti, Ø3.5 x 24mm TX, S/I Extrem	10
	15.14.03.35026	Locking Screw, Ti, Ø3.5 x 26mm TX, S/I Extrem	10
	15.15.03.27008	Headless Screw, Ti, Ø2.7 x 8mm TX, S/I Extrem	5
	15.15.03.27010	Headless Screw, Ti, Ø2.7 x 10mm TX, S/I Extrem	5
	15.15.03.27012	Headless Screw, Ti, Ø2.7 x 12mm TX, S/I Extrem	5
	15.15.03.27014	Headless Screw, Ti, Ø2.7 x 14mm TX, S/I Extrem	5
	15.15.03.27016	Headless Screw, Ti, Ø2.7 x 16mm TX, S/I Extrem	5
	15.15.03.27018	Headless Screw, Ti, Ø2.7 x 18mm TX, S/I Extrem	5
	15.15.03.27020	Headless Screw, Ti, Ø2.7 x 20mm TX, S/I Extrem	5
	15.15.03.35008	Headless Screw, Ti, Ø3.5 x 8mm TX, S/I Extrem	5
	15.15.03.35010	Headless Screw, Ti, Ø3.5 x 10mm TX, S/I Extrem	5
	15.15.03.35012	Headless Screw, Ti, Ø3.5 x 12mm TX, S/I Extrem	5
	15.15.03.35014	Headless Screw, Ti, Ø3.5 x 14mm TX, S/I Extrem	5
	15.15.03.35016	Headless Screw, Ti, Ø3.5 x 16mm TX, S/I Extrem	5
	15.15.03.35018	Headless Screw, Ti, Ø3.5 x 18mm TX, S/I Extrem	5
	15.15.03.35020	Headless Screw, Ti, Ø3.5 x 20mm TX, S/I Extrem	5








Implant Case Set 1













	N#	REF	DESC	QTY
	7	26.26.35.27352	Lock-Plate, Ti, Hallux Valgus, Sma, 0, IE 2,7/3.5	2
	8	26.26.35.27204	Lock-Plate, Ti, Hallux Valgus, L, Rig, 4, IE 2,7mm	2
	9	26.26.35.27104	Lock-Plate, Ti, Hallux Valgus, L, Lef, 4, IE 2,7mm	2
	10	26.26.35.27202	Lock-Plate, Ti, Hallux Valgus, L, Rig, 2, IE 2,7mm	2
	11	26.26.35.27102	Lock-Plate, Ti, Hallux Valgus, L, Lef, 2, IE 2,7mm	2
	12	26.26.35.27200	Lock-Plate, Ti, Hallux Valgus, L, Rig, 0, IE 2,7mm	2
	13	26.26.35.27100	Lock-Plate, Ti, Hallux Valgus, L, Lef, 0, IE 2,7mm	2
	14	09.05.15.00035	Drill Guide Impactor HV - 3,5 Straight	1
	15	26.26.35.35006	Lock-Plate, Ti, Hallux Valgus, 2P 3P, 6, IE 2,7/3.5	2
	16	26.26.35.35004	Lock-Plate, Ti, Hallux Valgus, 2P 3P, 4, IE 2,7/3.5	2
	17	26.26.35.35002	Lock-Plate, Ti, Hallux Valgus, 2P 3P, 2, IE 2,7/3.5	2
	18	26.26.35.27003	Lock-Plate, Ti, Hallux Valgus, Medium, IE 2,7mm	2




















Implant Case
Set 2

	N#	REF	DESC	QTY
	19	09.05.15.00000	Drill Guide/Impactor, HV Plate, Small	1
	20	26.26.35.27002	Lock-Plate, Ti, Hallux Valgus, Small, IE 2,7mm	2
	21	26.26.35.27001	Lock-Plate, Ti, Hallux Valgus, XSmall, 2, IE 2,7mm	2
	22	26.26.35.27000	Lock-Plate, Ti, Hallux Valgus, XSmall, 0, IE 2,7mm	2
	23	09.05.15.00040	Drill Guide/Impactor, HV Plate, 2.7/3.5 Curved	1
	24	26.26.35.27004	Lock-Plate, Ti, Hallux Valgus, Large, IE 2,7mm	1
	25	26.26.35.27357	Lock-Plate, Ti, Hallux Valgus, Lar, 7, IE 2,7/3.5	1
	26	26.26.35.27355	Lock-Plate, Ti, Hallux Valgus, Lar, 5, IE 2,7/3.5	1
	27	26.26.35.27353	Lock-Plate, Ti, Hallux Valgus, Lar, 3, IE 2,7/3.5	1
	28	26.26.35.27350	Lock-Plate, Ti, Hallux Valgus, Lar, 0, IE 2,7/3.5	1



	N#	REF	DESC	QTY
	29	09.01.07.20021	Drill Bit \varnothing 2.0 x 125mm, S-50, Stryker, B-4.5, Blue	2
	30	09.01.07.25021	Drill Bit \varnothing 2.5 x 125mm, S-50, Stryker, B-4.5, Yellow	2
	31	09.01.03.20020	Drill Bit \varnothing 2.0 x 120mm, S-50, AO, Blue	2
	32	09.01.03.25020	Drill Bit \varnothing 2.5 x 120mm, S-50, AO, Yellow	2
	33	09.02.01.27010	Countersink, \varnothing 2.7mm, Hallux Valgus, AO	2
	34	09.01.20.20015	Drill Bit 2in1, \varnothing 2.0mm, Stop 15mm, 90mm, Blue	2
	34	09.01.20.25015	Drill Bit 2in1, \varnothing 2.5mm, Stop 15mm, 90mm, Yellow	2
	35	09.01.17.20150	Graduated Drill Bit \varnothing 2,0 x 150 mm, AO Coupling, Blue	2
	35	09.01.17.25150	Graduated Drill Bit, \varnothing 2.5 X 150mm, AO coupling, Yellow	2

Box Instruments

	N#	REF	DESC	QTY
	36	09.34.01.00035	Preparation Chisel, HV Plate, 3.5	1
	37	09.34.01.00090	Preparation Chisel, Halux Valgus Plate, 90°	1
	38	09.34.01.02735	Preparation Chisel, HV Plate, Curved 2.7/3.5	1
	39	09.34.00.20004	Chisel 200mm, 4x3mm	1
	40	09.10.06.00150	Plate and Screw Holding Forceps, Angled, 150mm	1
	41	09.33.00.27220	HAJEK Hammer, 140g, Ø27 x 220mm	1
	42	09.08.01.00041	Depth Gauge 40mm, Long	1
	43	09.07.04.08091	Shaft Screwdriver, TX-8, 90mm, AO, Blue	2
	44	09.07.04.10092	Shaft Screwdriver, TX-10, 90mm, AO, Yellow	2
	45	09.06.00.46050	Drill Sleeve Ø4.6 x 50mm	2
	46	09.05.14.03020	Graduated Drill Guide, Threaded, Ø2,0x 30mm, Blue	1
	47	09.05.14.03025	Graduated Drill Guide, Threaded, Ø2,5x 30mm, Yellow	1
	48	09.20.03.00000	Fixation Wrench, Hallux Valgus Plate	1
	49	09.20.03.00035	Fixation Wrench, Hallux Valgus Plate 3.5mm	1
	50	09.04.04.12030	Cannulated Handle, AO, 120mm, Green	1
	51	09.05.15.02735	Drill Guide/Impactor, HV Plate, Small 2.7/3.5	1
	52	09.05.15.00010	Drill Guide/Impactor, HV Plate, Medium-Large	1
	53	09.05.15.00020	Drill Guide Impactor HV - 2,7 L	1





AST.BOX.025
REV.02
11/2022